

PRODUCT DETAIL



Ph 0800 27 26 26

info@ansutek.co.nz
www.ansutek.co.nz
www.ansutekbiz.co.nz

Wrapped connector

Can be used with various ranges of plastic and glass pipettes from 0.1 to 100ml (with or without graduations, including Pasteur pipettes)



More Effective

Comfortable grip

Widely used in clinical diagnostic experiments, biological experiments, pharmacy and chemical experiments, etc.



Easily Data Processing

Two drainage modes

The mode can be freely selected during the draining process:
[gravity draining mode](#) & [electric draining mode](#)

(Note: Regardless of the method, the speed can be adjusted sensitively by finger pressure)



Further pressing the drain button will trigger the micro pump to speed up the drain.

Humanization Design

Membrane protection

Adapt to 0.45µm hydrophobic filter membrane, hydrophobic high-temperature disinfection filter membrane to prevent damage to the equipment due to suction.



Humanization Design

Charging is more convenient

Strong resistance to over-discharge and high-speed discharge; safe to use, high power density



Charge for 2-3 hours, could work intermittently for about 8 hours